




Technical specification:

- Six-channel head, made of aluminium and powder coated in RAL 9002
- Front, rear cover powder coated in any color (RAL color palette)
- Head height: 600 mm
- Head depth: 374 mm
- Head width: 244 mm

Equipment options:

- Medical gas and vacuum outlets code types: AGA / DIN / BS / AFNOR
- Electrical sockets standards: Schuko / Legrand / BS / others
- IT sockets: RJ-11 / RJ-12 / RJ-45 of cat.5e or higher
- Nurse call systems of different kinds / producers
- LED lighting on the arm or beneath the head
- Pneumatic / electromagnetic brakes locking arms rotation

Sample equipment of Anesthesia pendant:

1. Medical gas outlets: 1x O2, 1x N2O, 2x AIR, 2x VAC, 1x AGSS
2. IP 44 electrical sockets: 8 pcs (230V/50Hz)
3. Equipotential sockets: 4 pcs
4. RJ45 sockets of cat. 6: 2 pcs
5. Rail (DIN 25x10mm type) – 300 mm
6. Arms pneumatic brake



Product characteristics

LISSAclassic-lift

Ceiling pendant
Version: 20170124ENG

Contact:

magda.baszak@inmed.pl
+71 398 02 28

www.inmed.pl



INMED-Karczewscy
LLC LLP
32a Kwiatowa Street
55-330 Krępice, Poland