



#### Technical specification:

- Six-channel head, made of aluminium and powder coated in RAL 9002
- Front, rear, top and bottom cover powder coated in any color (RAL color palette)
- Two vertical supporting tubes from the front
- Head height: 600/700/800/900/1000/1100/1200 mm
- Head depth (with supporting tubes): 348 mm
- Head width: 270 mm

#### Equipment options:

- Medical gas and vacuum outlets code types: AGA / DIN / BS / AFNOR
- Electrical sockets standards: Schuko / Legrand / BS / others
- IT sockets: RJ-11 / RJ-12 / RJ-45 of cat.5e or higher
- Nurse call systems of different kinds / producers
- Shelves, medical rails, IV poles, catheters baskets or wound care products carriers
- LED lighting on the arm or beneath the head
- Pneumatic / electromagnetic brakes locking arms rotation

#### Sample equipment of Surgical pendant:

1. Medical gas outlets: 2x AIR, 2x VAC, 1x AIR-8, 1x CO2
2. IP 44 electrical sockets: 12 pcs (230V/50Hz)
3. Equipotential sockets: 12 pcs
4. RJ45 sockets of cat. 6: 2 pcs
5. Provisions for User's Installation: 3 pcs
6. IV pole: 1pc
7. Shelf with side rails (DIN 25x10mm type): 1 pc
8. Pneumatic or electromagnetic arm brake



#### Product characteristics

### LISSAclassic

Ceiling pendant  
Version: 20170918ENG

#### Contact:

magda.baszak@inmed.pl

+71 398 02 28

www.inmed.pl



INMED-Karczewscy  
LLC LLP  
32a Kwiatowa Street  
55-330 Krępiec, Poland