


hidden expectations

MERY cabinetto






Hidden expectations

MERYcabinetto is the only unit on the market accommodating bed head unit along with an additional bed for the accompanying person, e.g. when the mother needs to stay in hospital with the child. The Cabinetto® system allows pulling out the matters with support from the cabinet that remains fully functional. The bed is pulled out automatically, in 30 seconds by using the button on the unit or remote control.

MERYcabinetto accommodates all necessary fittings as medical gas outlets, electrical HV and LV services along with lighting. The fittings can be located on a dedicated panel (mBox), on the side of the cabinet or on the wall depending on the customers' requirements. Each unit can be customized to match the room's interior. MERYcabinetto is built from the highest quality CORIAN® by DuPont® boards with the wide range of finishing components. If desired it can be also built from non-ignitable (Class B) laminated boards. Wood grain can be customized.

The cabinets come in a variety of shapes and sizes and accommodate different functionalities as self-closing door hinges and drawers, coat hangers or case for bedding. Size of the mattress 80x200cm.





Design flexibility

MERYcabinetto can meet any demands as far as the style of the patients' room is concerned. Moreover, it can easily transform any hospital's room interior, especially the rooms for the little ones. We provide CORIAN® by DuPont® colourful panels, a variety of patterns or graphics to build fancy units for kids' rooms and wards.

And because we know how important for a kid is to have a parent around while staying in a hospital – we have given the parents the opportunity to stay in the same room with the kid and rest on comfortable full-size bed.

