



#### Specification:

- Aluminium frame mounted directly on the wall
- mBOX - consisting of 3 independent ducts to plug in all necessary media
- Horizontal lighting frame made of brushed aluminium
- AmbientLight located on the frame
- Cabinet with folding bed for accompanying person
- Length: single unit with cabinet – 2430mm (1500 + 930mm)
- Height of cabinet: 2480mm
- Depth with cabinet: 680 mm: depth of unit above bed – 128 mm
- Dimensions with folded bed: ~2600mm

#### Equipment options:

- Medical gases and vacuum outlets code types: AGA / DIN / BS / AFNOR
- Lighting: T5 / KLD / LED
- Electric sockets standards: German / French / British / others
- IT sockets: RJ-11 / RJ-12 / RJ-45 of cat.5e or higher
- Nurse call system
- Lighting switches on patient's handset
- Decorative covers material: Corian, laminated particle board



#### Standard equipment for hospital rooms:

1. Medical gases outlets: 2 x (O2 and VAC)
2. Electrical sockets: 3 x (230V/50Hz)
3. Equipotential socket: 1 x
4. RJ45 cat.6 socket for IT: 1 x
5. Mounting holes for nurse call system
6. Indirect light 2x28W: switch on the unit
7. Direct light 1x28W: switch on the unit
8. LED AmbientLight 2x4,8W LED: switch on the unit

#### Product characteristics

## MERYcabinetto

Bed head units

Version: 20151127ENG

Contact: Magda Baszak  
 magda.baszak@inmed.pl  
 www.inmed.pl



INMED-Karczewscy  
 LLC LLP  
 Kwiatowa 32A  
 55-330 Kępice