



#### Specification:

- Body made of anodized aluminium profile mounted directly on the wall
- Front cover powder coated in any color (RAL color palette)
- 1 electrical and 1 gas installation input for double and multi-units
- Length: single unit-2000 mm, double unit-4000 mm, triple unit-6000 mm
- Depth: 110 mm
- Height: 235 mm

#### Equipment options:

- Body color: anodized / powder coated in any color (RAL color palette)
- Three, four – chamber body
- Medical gases and vacuum outlets code types: AGA / DIN / BS / AFNOR
- Lighting: T5 / KLD / LED
- Electric sockets standards: German / French / British / other
- IT sockets: RJ-11 / RJ-12 / RJ-45 of cat.5e or higher
- Nurse call system
- All components – subject to being touched by hospital staff in managing the patient may be manufactured from antimicrobial copper (AMC)
- Additional equipment: garbage for disposable items, drip lamps, shelves of various dimensions and carrying capacity, dual-swivel monitor arms, infusion liquids jibs extension arms or medical mobile trolleys with storage shelves and drawers

#### Standard equipment for Intensive Care Units:

1. Medical gases outlets: 6 x (2 x O<sub>2</sub>, 2 x AIR, 2 x VAC)
2. Electrical sockets: 8 x (230V/50Hz)
3. Equipotential sockets: 4 x
4. RJ45 cat.6 sockets for IT: 2 x
5. Mounting holes for nurse all system
6. Indirect light 2x54W T5: switch on the unit
7. Direct light 1x36W KLD: switch on the unit
8. Night light 1x2,8W LED: switch on the unit

#### Medical rails – options:

1. Mounted on the unit – 1 pc
2. Mounted on the unit – 2 pcs
3. Mounted on the wall – 1 pc
4. Mounted on the wall – 2 pcs

#### Product characteristics

**MERYcare**

Bed head unit

Version: 20170123ENG

Contact: Magda Baszak  
magda.baszak@inmed.pl  
www.inmed.pl



INMED-Karczewscy  
LLC LLP  
Kwiatowa 32A  
55-330 Krepice