



Specification:

- Body made of anodized aluminium profile mounted directly on the wall
- Front cover powder coated in any color (RAL color palette)
- Length: single unit-2000 mm, double unit-4000 mm, triple unit-6000 mm
- Depth: 65 mm
- Height: 300 mm
- 1 electrical and 1 gas installation input for double and multi-units

Equipment options:

- Body color anodized / powder coated in any color (RAL color palette)
- Four, five – chamber body
- Medical gases and vacuum outlets code types: AGA / DIN / BS / AFNOR
- Lighting: T5 / KLD / LED
- Electric sockets standards: German / French / British /other
- IT sockets: RJ-11 / RJ-12 / RJ-45 of cat.5e or higher
- Nurse call system
- All components – subject to being touched by hospital staff in managing the patient may be manufactured from antimicrobial copper (AMC)
- Additional equipment: garbage for disposable items, drip lamps, shelves of various dimensions and carrying capacity, dual-swivel monitor arms, infusion liquids jibs extension arms or medical mobile trolleys with storage shelves and drawers



Standard equipment for Intensive Care Units:

1. Medical gases outlets: 6 x (2 x O2, 2 x AIR, 2 x VAC)
2. Electrical sockets: 6 x (230V/50Hz)
3. Equipotential sockets: 2 x
4. RJ45 cat.6 sockets for IT: 2 x
5. Mounting holes for nurse call system
6. Indirect light 2x36W KLD: switch on the unit
7. Direct light 1x36W KLD: switch on the unit
8. Night light 1x2,8W LED: switch on the unit

Medical rails – options:

1. Mounted on the unit – 1 pc
2. Mounted on the unit – 2 pcs
3. Mounted on the wall – 1 pc
4. Mounted on the wall – 2 pcs

Product characteristics

MERYcare-light

Bed head unit

Version: : 20170123ENG

Contact: Magda Baszak
 magda.baszak@inmed.pl
 www.inmed.pl



INMED-Karczewscy
 LLC LLP
 Kwiatowa 32A
 55-330 Krępiec