


E – Power supply output  
G – Gas supply output

Entrance  
←



#### Specification:

- Body made of aluminium profile mounted directly on the wall
- Front decorative cover made of brushed or anodized aluminium
- Length: single unit 1500 mm
- Depth: 160 mm
- Height: 164 mm

#### Equipment options:

- Body color: anodized / powder coated in any color (RAL color palette)
- Medical gases and vacuum outlets code types: AGA / DIN / BS / AFNOR
- Lighting: T5 / LED
- Electric sockets standards: Schuko / Legrand / BS / others
- IT sockets: RJ-11 / RJ-12 / RJ-45, cat.5e or higher
- Nurse call system
- Lighting switches on patient's handset

#### Standard equipment for hospital rooms:

1. Medical gases outlets: 2 x (O<sub>2</sub> and VAC)
2. Electrical sockets: 3 x (230V/50Hz)
3. Equipotential socket: 1 x
4. RJ45 cat.6 socket for IT: 1 x
5. Mounting holes for nurse call system
6. Direct light 1x39W T5: switch on the unit
7. Direct side-light 1x39W T5: switch on the unit
8. LED Ambient Light 1x2,8W LED: switch on the unit

#### Product characteristics

**MERYedge**

Bed head unit  
Version: 20170123PL

#### Contact:

magda.baszak@inmed.pl  
+71 398 02 28

www.inmed.pl



INMED-Karczewscy  
Sp. z o.o. Sp. K.  
ul. Kwiatowa 32a  
55-330 Kępice