


**E** – Power supply output  
**G** – Gas supply output



**Specification:**

- Body made of anodized aluminium profile mounted directly on the wall
- Front cover powder coated in any color (RAL color palette)
- Length: single unit-1600 mm, double unit-3200 mm, triple unit-4800 mm
- Depth: 65 mm
- Height: 300 mm
- 1 electrical and 1 gas installation input for double and multi-units

**Equipment options:**

- Body color: anodized/ powder coated in any color (RAL color palette)
- Three, four or five-chamber body
- Medical gases and vacuum outlets code types: AGA / DIN / BS / AFNOR
- Lighting: T5 / KLD / LED
- Electric sockets standards: German / French / GB / others
- IT sockets: RJ-11 / RJ-12 / RJ-45 of cat.5e or higher
- Nurse call system
- Lighting switch on patient's handset

**Standard equipment for hospital rooms:**

1. Medical gases outlets: 2 x (O2 and VAC)
2. Electrical sockets: 3 x (230V/50Hz)
3. Equipotential socket: 1 x
4. RJ45 cat.6 socket for IT: 1 x
5. Mounting holes for nurse call system
6. Indirect light 2x36W KLD: switch on the unit
7. Direct light 1x24W KLD: switch on the unit
8. Night light 1x2,8W LED: switch on the unit

**Product characteristics**

**MERYlight**

Bed head unit

Version: 20170123ENG

Contact: Magda Baszak  
 magda.baszak@inmed.pl  
 www.inmed.pl



INMED-Karczewscy  
 LLC LLP  
 Kwiatowa 32A  
 55-330 Krepice