

wall-mounted

# MERYwall

DI: MERY20140325EN



## the inspiration of future

Our constant challenge has always been the search for inspiration to release new product. We all have been totally astonished with the result of the search for inspiration which is functionality and beauty blended together. The employment of wood, not commonly used in manufacture of medical products and aluminum, resulted in outrageous medical supply unit that, for a precious moment, makes hospitalized patients feel at home.

MERYwall supply units are in compliance with the philosophy of unique and outrageous solutions of in-hospital rooms' arrangement. We are proud of our variety of colors available and ergonomic functions combined together to contribute to modern and convenient design of in-hospital rooms.

MERYwall turns to be the result of long tradition and experience of INMED blended with visions inspired by future.





# the perfection of light

*„Never fear shadows, for shadows only mean there is a light shining somewhere near by.“*  
*Oscar Wilde*

Light activates our imagination, directing us to the areas we do not regularly notice. MERYwall is a unique combination of light symbolizing functionality and decorativeness. For us light turns to be much more than a component of technical equipment, for us light equals perfection and our focus is to make it excellently perfect.

Our product is based on esthetic, horizontal lighting frame made of brushed aluminum, perfectly suiting rest of equipment. Among others the source of light is line fluorescent light activated by electronic ballast. Modern design has resulted in unique visual effect. The light controlling options in the frame may be adjusted individually. Patients may also be equipped with manipulator to control the light. Unique design and exceptional LED functionality is represented by AmbientLight located in the frame. Thanks to it the patient may recover at nice, delicate and soothing light.



Full light



Direct light



AmbientLight










# colors

For us confronting colors with health and safety has turned out to be real challenge. MERYwall was manufactured of flame-resistant, laminated panel B class of fire or Corian® by DuPont®. The most sophisticated demands and styles of the customers are met thanks to the wide variety of textures and structured surfaces of the panels. Apart from fire-resistance and sophisticated design each panel is granted E1 hygienic certificate (in compliance with Polish standards). We particularly appreciate the fact that all products granted E1 class are human and environment friendly and recommended in health services.

If the design of the room does not include wooden elements, we offer Corian® by DuPont®, balancing unique beauty and functionality. It is particularly applicable to innovative and creative long-lasting designs. Perfectly smooth and clear, polished to a high gloss, with no bonds visible, commonly used in operating theaters has recently been highly recommended in in-hospital rooms.

Patterns and colors available:



# flexibility of solutions

Clear and easy electric energy supply and access to medical gas sources are outstanding features of MERYwall. mBOX, consisting of 3 independent ducts to plug in all necessary media, is located on a decorative wooden panel made of fire-resistant materials or Corian®. For higher quality media management the ducts have been separated into electrical components duct for plug-ins and light switching, various medical gas outlets, and IT duct. Each is equipped with easily accessible, accessories friendly and spacious area to fit own selection of equipment needed.

MERYwall does not only combine media in one location but also provides an easy access to all accessories required by patient or hospital staff. mBOX front panel is made of brushed aluminum, however to let the cat out of the bag a little, we have been recently focusing on new technologies to be applied in mBOX manufacture.





# adjustability

MERYwall means uniqueness of each solution. Every in-hospital room can have exceptional atmosphere by application of wide range of colors, equipment, shapes of mBoxes or light. The number of configurations is only limited by designer imagination. In fact every piece of additional equipment like elements of beds, chairs, furniture may correspond to philosophy of MERYwall together succeeding in perfect adjustability.

MERYwall is dedicated to single and multi-bed solutions ensuring high level of ergonomics and unique design. High quality of material employed in manufacture of MERYwall provides the long lasting highest standard of in-hospital rooms.

