

comfortable solutions


MERY comfort



Comfort and ergonomics

MERYcomfort is one of the few units on the market that was designed to be as much functional and accessible for the patients as possible. All the fittings the patients need – sockets, light switches, USB ports for charging mobile devices and nurse call buttons – are accommodated within a hand's reach.

The unit provides an innovative direct lighting in double or more beds versions. The light is located in a vertical part of the bed head and directed at the particular bed. Besides standard fittings (medical gas outlets, electrical sockets, USB sockets, lighting) MERYcomfort provides shelves that are integrated with the body. The shelves' load capacity depends on required devices.









Components

It is always highly desirable that a bed head matches the interior design, mostly as far as the colour is concerned. MERYcomfort goes much further. The style, colour and materials can be customized. The following materials can be used in manufacturing:

- CORIAN® by DuPont®
- Non-ignitable (Class B) laminated boards (wood grain can be customized)
- Toughened glass with sand-blasted or printed décor
- Plastic or other safe materials with printed décors

Individuality

Our main focus has always been to organize the space for patients as well as possible and make it as personalized as possible. MERYcomfort is built from aluminium frame that enables transferring the loads and fixing different equipment. INMED takes pride in supplying diverse MERYcomfort units developed with cooperation and involvement of our customers, e.g. the unit with integrated cooler for medicines, the unit with neon decorative lighting or with the integrated shelf. We treat every project as a challenge considering the unlimited possibility of arrangement. Each MERYcomfort is a big challenge for us as we face with an unlimited number of possible arrangements. The vast majority of these units are absolutely unique. Our engineers together with the customers work on every project until they achieve the desired effect. We are awaiting the new ideas and are ready to face new challenges. We are certain the final effect will be astonishing. As always.

see a movie about this product



Boxes integrated with the unit



Shelves for medical instruments and holders for patient's mobile devices



