


elegance in simplicity

MERYedge





'Fashion changes, but style endures'
Gabrielle Bonheur Chanel

Inspired by vicinity

MERYedge is the innovative solution for patients' rooms. This currently is the only bed head unit on the market that is constructed as a one single block and do not require assembling and disassembling part by part. MERYedge style is characterized by a unique simplicity and elegance. The rooms this unit is installed in gain a completely new character and an unconventional, unusual look.

Construction solutions applied enable hiding all electrical and LV sockets along with lampshades what gives the impression of simplicity, transparency and timelessness.

MERYedge is not only a bed head unit but also the element of bigger arrangements.








Engineering ingenuity

MERYedge accommodates all services required – medical gas outlets, electrical HV and LV sockets and different types of lighting. The design, construction and ergonomic arrangement of electrical sockets at the bottom of the unit result in its exceptional style.

However first of all MERYedge has the following advantages:

- Durable, anodized aluminium 2nd fix of gas outlets
- Brushed aluminium body
- One-piece construction
- Dedicated duct for accessories as shelves, trays, holders etc.



movable accessories



direct light and room light